


# QUART didoMAS

## Advanced mA/mAs and Time Measurements



# QUART didoMAS

Advanced Technology and Improved Methodology



The QUART didoMAS meter is designed for easy and very precise current and current-time product measurements.

The QUART didoMAS meters automatically set the range of measurement. No pre-setting is required for direct reading of mA, mAs and time parameters.

The meters can be used throughout the complete range of radiographic equipment including fluoroscopic or mammography exposures.

The method of measurement is invasive into the circuit. That delivers the most reliable and reproducible results.

The didoMAS is equipped with a backlit display to assure swift readings even in darkened environments. To provide the ability to track generator characteristics, the mAs is refreshed continuously on the meter display while the measurement is running.

The meter is powered by a rechargeable battery. One charge is sufficient to last approximately 80 hours of continuous use. Recharging the meter until full takes only between 3 - 4 hours. A warning will appear on the display when the battery is running low.

## Special Features

The meter features an extra-long cable between the base and the detector unit. A customised 1.5 m long cable for the connection between detector head and circuit is included in the delivery.


The connection between circuit and the detector unit is polarity independent. The mA is refreshed and displayed four times per second.

## TECHNICAL SPECIFICATIONS

Temp. Range	15 - 30°C	
Storage Range	-10 - 65°C	
Environmental	Humidity 20-75% / Air Pressure 20g/m <sup>3</sup>	
Weight	Base unit:	180 g
	Detector unit:	120 g including cable
Size	Base unit:	17.0 x 7.0 x 4.5 cm (LxWxH)
	Detector unit:	4.7 x 3.0 x 4.0 cm (LxWxH)

## QUART didoMAS

mAs	Range:	0.01 - 9999 mAs
	Uncertainty:	± 1%
mA	Range:	0.1 - 800 mA
	Modes (3):	<ul style="list-style-type: none"> <li>• One-second-mean current refreshed four times per second (real-time)</li> <li>• Maximum mean current (&gt; 1s)</li> <li>• Total average current</li> </ul>
	Uncertainty:	± 1%
	Exposure Time	Range:
	Mode:	Duration of current flow
	Uncertainty:	± 1% (for rectangular signal above 50 ms)



 **Active Radsys**  
[www.activeradsys.it](http://www.activeradsys.it) [info@activeradsys.it](mailto:info@activeradsys.it)  
 Tel: 0544 408071 Fax: 0544 201477

 **QUART**  
 Quality Assurance in Radiological Technologies