






QUART DVT_TEC Automatic Image Evaluation Software

Accompanies QUART CBCT/DVT_AP phantom (complies DIN 6868-161)



- Unique solution for CBCT / Digital Volume Tomography QA/QC procedures
- Auto-calculation of Nyquist Frequency, Homogeneity and Contrast-to-Noise Ratio
- Auto-calculation of Modulation Transfer Function
- Auto-calculation of System Indicator to provide system inter-comparability
- Relate imaging performance to equipment specific dose parameters (methodology included)
- Time-saving, comprehensive imaging performance assessment for CBCT/DVT
- Save data for future reference or archive purposes
- Print customer-specific protocols/hardcopy to document test results